

# Peptide Reagents



The creation of a peptide linkage between two amino acid segments is one of the most important reactions in organic and bioorganic chemistry. Many different strategies have been devised for selective amide bond formation from a carboxylic acid and an amino group, usually involving protection, activation, coupling and deprotection steps. The reagents listed in this section have been used in one or other of these steps in the synthesis of peptide molecules. An enormous body of literature exists on the subject of peptide synthesis, of which a selection of books and review articles is listed.<sup>1-12</sup>

## Protecting group reagents

To ensure specific coupling between the required carboxyl and amino groups, a range of protecting groups have been developed which can be selectively introduced and removed. References to methods of protection and deprotection can be found under specific items in the main Catalogue. Further information on the use of protecting groups in peptide synthesis can be found in the general references on peptide synthesis and in other specialist sources.<sup>13-15</sup>

## CARBOXYL PROTECTION

Although a wide variety of esters and other groups have been used for the protection of carboxylic acids, only a limited number of these find significant use in peptide synthesis.

### tert-Butyl ester

**L08855** tert-Butyl acetate

**L00431** N,N-Dimethylformamide di-tert-butyl acetal

**B22039** tert-Butyl 2,2,2-trichloroacetimidate

**L12338** tert-Butanol

### 1-Adamantyl ester

**A10209** 1-Adamantanol

### Dicyclopropylmethyl ester

**L01816** Dicyclopropylmethanol

### Benzylic esters

**A15188** 9-Anthracenemethanol

**L03292** Benzyl alcohol

**A13535** Benzyl bromide

**A12481** Benzyl chloride

**A13579** 9-(Chloromethyl)anthracene

**A12859** 4-(Chloromethyl)pyridine hydrochloride

**A15212** 9-Fluorenylmethanol

**A15559** 4-Methoxybenzyl alcohol

**A15742** 4-Nitrobenzyl alcohol

**A13127** 2-Nitrobenzyl bromide

**A15236** 4-Nitrobenzyl bromide

**A15749** 4-Nitrobenzyl chloride

**L00555** 2,3,4,5,6-Pentamethylbenzyl chloride

**A14698** Pyridine-4-methanol

**L00423** 2,4,6-Trimethylbenzyl chloride

## Alfa Aesar Worldwide Sales Offices

### NORTH AMERICA

Tel: 1-800-343-0660 or  
1-978-521-6300

Fax: 1-800-322-6350  
Email: info@alfa.com

### GERMANY

Tel: 00800 4566 4566 or  
+49-721-84007-280

Fax: 00800 4577 4577 or  
Email: Eurosales@alfa.com

### UNITED KINGDOM

Tel: 0800-801812 or  
+44 (0)1524-850506

Fax: +44 (0)1524-850608  
Email: UKsales@alfa.com

### FRANCE

Tel: 0800 03 51 47 or  
+33 03 88 62 26 90

Fax: +33 03 88 62 26 81  
Email: frventes@alfa.com

### INDIA

Tel: +91 (0)44 3251 9881  
Fax: +91 (0)44 2814 4154

Email: IndiaSales@alfa.com

### CHINA

Tel: (010) 51660581  
Fax: (010) 51987480

Email: sales-china@alfa.com

# Peptide Reagents Con't

## Benzhydryl-type ester precursors

**A12266** Benzophenone hydrazone

**B21706** 9-Fluorenone hydrazone

## Phenacyl esters

**A15576** 2-Bromoacetophenone

**A13415** 2-Bromo-4'-methoxyacetophenone

## Allyl ester

**A15026** Allyl alcohol

**A11766** Allyl bromide

## 2-Substituted ethyl esters

**A10275** 2-Bromoethanol

**L19365** 2-(p-Toluenesulfonyl)ethanol

**B23579** 2-(2-Hydroxyethyl)pyridine

**L08163** 2,2,2-Trichloroethanol

**L05442** 2-(Methylthio)ethanol

**B20970** 2-(Trimethylsilyl)ethanol

## Miscellaneous carboxyl protecting reagents

**A13005** Benzyl carbazate

**A12018** N-(Chloromethyl)phthalimide

**A12383** tert-Butyl carbazate

**B21292** N-(Hydroxymethyl)phthalimide

**A15238** 2-Chloroacetamide

## AMINO PROTECTION

The most popular protecting groups for the amino function are carbamates, particularly Benzyloxycarbonyl [Cbz, Z], Boc and Fmoc groups, but many other groups including trifluoroacetyl and trityl can be used.

## Benzyl carbamate [Cbz, Z]

**A15682** Benzyl chloroformate

**A12153** N-(Benzyloxycarbonyloxy)succinimide

## Substituted benzyl carbamates

**B25632** 4-Nitrobenzyl chloroformate

## tert-Butyl carbamate [Boc]

**A14708** Di-tert-butyl dicarbonate

**A18906** Boc-ON [2-(tert-Butoxycarbonyloximino)]

**L00506** N-Boc-imidazole 2-phenylacetonitrile]

**B22435** S-Boc-2-mercapto-4,6-dimethylpyrimidine

**A12383** tert-Butyl carbazate

## Benzo[b]thiophenesulfoylmethyl carbamate [Bsmoc]

**L19559** 1,1-Dioxobenzo[b]thiophen-2-ylmethyl chloroformate

**L19733** 1,1-Dioxobenzo[b]thiophen-2-ylmethyl N-succinidyl carbonate

## 2,2,2-Trichloroethyl carbamate [Troc]

**L06875** 2,2,2-Trichloroethyl chloroformate

### Alfa Aesar Worldwide Sales Offices

#### NORTH AMERICA

Tel: 1-800-343-0660 or  
1-978-521-6300  
Fax: 1-800-322-6350  
Email: info@alfa.com

#### GERMANY

Tel: 00800 4566 4566 or  
+49-721-84007-280  
Fax: 00800 4577 4577 or  
Email: Eurosales@alfa.com

#### UNITED KINGDOM

Tel: 0800-801812 or  
+44 (0)1524-850506  
Fax: +44 (0)1524-850608  
Email: UKsales@alfa.com

#### FRANCE

Tel: 0800 03 51 47 or  
+33 03 88 62 26 90  
Fax: +33 03 88 62 26 81  
Email: frventes@alfa.com

#### INDIA

Tel: +91 (0)44 3251 9881  
Fax: +91 (0)44 2814 4154  
Email: IndiaSales@alfa.com

#### CHINA

Tel: (010) 51660581  
Fax: (010) 51987480  
Email: sales-china@alfa.com

# Peptide Reagents Con't



## 9-Fluorenylmethyl carbamate [Fmoc]

**A11683** 9-Fluorenylmethyl chloroformate

**A13143** N-(9-Fluorenylmethoxycarbonyloxy)-succinimide

## Allyl carbamate [Alloc]

**B22332** Allyl chloroformate

## Trifluoroacetyl group

**A13614** Trifluoroacetic anhydride

**L00135** Trifluoroacetylimidazole

## Trityl groups (Trt, Mmt)

**A12864** Bromotriphenylmethane

**A11799** Chlorotriphenylmethane

**A10545** 4-Methoxytrityl chloride

## Miscellaneous amino protecting groups

**A14822** N-(Ethoxycarbonyl)phthalimide

**L06432** Chlorocarbonylsulfonyl chloride

**L00770** 2-Chloro-3,5-dinitropyridine

**B23627** 2-Nitrobenzenesulfonyl chloride

**A15857** 3-Nitrophthalic anhydride

## SIDE CHAIN PROTECTION

The presence of reactive side chains in certain amino acids can interfere with peptide synthesis. Some examples are given here of reagents used to mask this reactivity:

### Hydroxyl protection

Of the vast range of groups available for the blocking of hydroxyl groups, only a few find regular use in peptide synthesis for serine, threonine and tyrosine derivatives. Benzyl ethers (from benzyl halides or alcohol) and tert-butyl ethers are widely used. For silyl protection, see Appendix 4. Among the most useful OH protecting reagents are:

**L13471** N-(2-Bromobenzyloxycarbonyloxy)-succinimide

**A12387** Benzyl 2,2,2-trichloroacetimidate

**A10370** 3-Bromocyclohexene

**L13576** 3-Bromopentane

**B22039** tert-Butyl 2,2,2-trichloroacetimidate

**A12859** 4-(Chloromethyl)pyridine hydrochloride

**L01050** 2-Methoxyethoxymethyl chloride

### Thiol protection

A variety of methods have been employed to mask the thiol group of cysteine, including benzyl, substituted benzyl, benzhydryl, trityl, Cbz and Boc groups.<sup>16</sup> Examples of reagents for SH protection are given:

**A12884** Benzhydryl

**A13820** N-(Hydroxymethyl)acetamide

**L14367** trans- $\beta$ -Nitrostyrene

**A10366** Triphenylmethanol

**L13316** 4-Vinylpyridine

## Alfa Aesar Worldwide Sales Offices

### NORTH AMERICA

Tel: 1-800-343-0660 or  
1-978-521-6300  
Fax: 1-800-322-6350  
Email: info@alfa.com

### GERMANY

Tel: 00800 4566 4566 or  
+49-721-84007-280  
Fax: 00800 4577 4577 or  
Email: Eurosales@alfa.com

### UNITED KINGDOM

Tel: 0800-801812 or  
+44 (0)1524-850506  
Fax: +44 (0)1524-850608  
Email: UKsales@alfa.com

### FRANCE

Tel: 0800 03 51 47 or  
+33 03 88 62 26 90  
Fax: +33 03 88 62 26 81  
Email: frventes@alfa.com

### INDIA

Tel: +91 (0)44 3251 9881  
Fax: +91 (0)44 2814 4154  
Email: IndiaSales@alfa.com

### CHINA

Tel: (010) 51660581  
Fax: (010) 51987480  
Email: sales-china@alfa.com

# Peptide Reagents Con't

## Carboxyl protection

The side chain carboxyl groups of aspartyl and glutamyl residues can be blocked by means of one of the reagents listed above for carboxyl protection, and also:

**A10697** 2-Adamantanol

**A10381** Pyridine-3-methanol

## Guanidino protection

The guanidino group of arginine has been protected by N-nitration, or by Cbz or Boc derivatives, amongst others. Arylsulfonyl protection has been found to be particularly useful:

**A11775** Mesitylenesulfonyl chloride

**L19561** 2,2,4,6,7-Pentamethyldihydro-

**L11829** 4-Methoxy-2,3,6-trimethylbenzenesulfonyl chloride

benzo[b]furan-5-sulfonyl chloride

## Imidazole protection

The imidazole ring in histidine can interfere with acylation reactions and can also promote racemization. It has been blocked by a variety of methods including formation of the Boc or Troc derivatives (reagents listed under Amino Protection) and the following:

**A11871** 2,4-Dinitro-1-fluorobenzene

**A13127** 2-Nitrobenzyl bromide


**A13788** 4-Methoxybenzenesulfonyl chloride

**A14547** p-Toluenesulfonyl chloride

## Activation and coupling methods

Formation of the peptide link under sufficiently mild conditions normally requires activation of the carboxylic acid function by conversion to a more electrophilic species, such as an acyl halide, azide, anhydride, mixed anhydride or active ester, which then undergoes coupling *in situ*, or as a separate step, with the amino function of the second component, as shown in Scheme 1.

### Scheme 1



Acid chlorides have limited value in peptide coupling because of the danger of hydrolysis, racemization, cleavage of protecting groups and other side reactions.<sup>12</sup> These difficulties can generally be avoided by the use of acid fluorides as active intermediates.

### Alfa Aesar Worldwide Sales Offices

#### NORTH AMERICA

Tel: 1-800-343-0660 or  
1-978-521-6300  
Fax: 1-800-322-6350  
Email: info@alfa.com

#### GERMANY

Tel: 00800 4566 4566 or  
+49-721-84007-280  
Fax: 00800 4577 4577 or  
Email: Eurosales@alfa.com

#### UNITED KINGDOM

Tel: 0800-801812 or  
+44 (0)1524-850506  
Fax: +44 (0)1524-850608  
Email: UKsales@alfa.com

#### FRANCE

Tel: 0800 03 51 47 or  
+33 03 88 62 26 90  
Fax: +33 03 88 62 26 81  
Email: frventes@alfa.com

#### INDIA

Tel: +91 (0)44 3251 9881  
Fax: +91 (0)44 2814 4154  
Email: IndiaSales@alfa.com

#### CHINA

Tel: (010) 51660581  
Fax: (010) 51987480  
Email: sales-china@alfa.com


# Peptide Reagents Con't

Mixed anhydrides of various kinds, formed by reaction of a carboxylic acid with chloroformates, acid chlorides, etc., in the presence of a tertiary amine base, have a long history as peptide intermediates. A variety of phosphorus-based reagents have been employed to form mixed anhydrides and some have shown promising results, often with simplified procedures.

However, for more than four decades the most widely-used reagents have been carbodiimides, especially N,N-dicyclohexylcarbodiimide [DCC] and water-soluble carbodiimide<sup>17</sup> [EDCI, "WSC"]. These will couple protected amino acids directly, but better results are often obtained with various additives or activating agents, most of which can form an "active ester", containing a good leaving group, with the carboxyl function. Frequently the coupling is performed as two separate steps: carbodiimide promoted formation (and possible isolation) of the active ester, followed by reaction with the free amino species.

The main problem is often partial loss of chirality caused by side reactions, the most important of which is considered to be azlactone (oxazolinone) formation by intramolecular cyclization of an N-acylated activated acid derivative with proton abstraction at the chiral center (Scheme 2).

**Scheme 2**



The suppression of racemization has been a major goal of much of the effort in development of coupling methods, and many techniques have been introduced, of which the use of active esters plays a prominent role. Newer coupling methods include the use of pyridinium salts, phosphonium salts and a variety of uronium salt reagents, some of which have been used to form active esters, giving superior results in specific coupling applications. Several other methods of peptide coupling are sometimes used, including the oldest of all, the azide method<sup>11</sup> due to Curtius,<sup>18</sup> of which the use of diphenylphosphonic azide<sup>19</sup> is a variation.

## Alfa Aesar Worldwide Sales Offices

### NORTH AMERICA

Tel: 1-800-343-0660 or  
1-978-521-6300  
Fax: 1-800-322-6350  
Email: info@alfa.com

### GERMANY

Tel: 00800 4566 4566 or  
+49-721-84007-280  
Fax: 00800 4577 4577 or  
Email: Eurosales@alfa.com

### UNITED KINGDOM

Tel: 0800-801812 or  
+44 (0)1524-850506  
Fax: +44 (0)1524-850608  
Email: UKsales@alfa.com

### FRANCE

Tel: 0800 03 51 47 or  
+33 03 88 62 26 90  
Fax: +33 03 88 62 26 81  
Email: frventes@alfa.com

### INDIA

Tel: +91 (0)44 3251 9881  
Fax: +91 (0)44 2814 4154  
Email: IndiaSales@alfa.com

### CHINA

Tel: (010) 51660581  
Fax: (010) 51987480  
Email: sales-china@alfa.com

# Peptide Reagents Con't



## ACTIVATING AND COUPLING REAGENTS

Further information and literature references on the use of most the following reagents can be found in the text entries in the main section of the Catalogue.

### Acyl halides

**L18013** Bis(tetramethylene)fluoroformamidinium hexafluorophosphate [BTFFH]

**A15666** Cyanuric fluoride

**A11992** Diethylaminosulfur trifluoride

**A18012** Oxalyl chloride

### Carbonic and carboxylic mixed anhydrides

**A14688** N,N'-Carbonyldiimidazole

**A13724** 2-Ethoxy-1-ethoxycarbonyl-1,2-dihydroquinoline [EEDQ]

**A14692** Isobutyl chloroformate

**L14159** 2,4,6-Trichlorobenzoyl chloride

**A15051** Trimethylacetyl chloride

### Phosphorus mixed anhydrides

**L08775** Bis(2-oxo-3-oxazolidinyl)phosphinic chloride [BOP-Cl]

**L09919** Diethyl chlorophosphite

**A11721** Diphenylphosphinic chloride

**A13546** Diphenyl phosphorochloridate

**A12724** Ethylene chlorophosphite

**A14530** Lawesson's Reagent

**L19271** 1-Propylphosphonic acid cyclic anhydride, [®T3P], 50+% soln. in DMF

**L11911** 50+% soln. in ethyl acetate

### Sulfonic mixed anhydride

**L12147** 3,5-Dichloro-2-hydroxybenzenesulfonyl chloride

### Carbodiimides

**L19463** N-Cyclohexylcarbodiimide on Merrifield resin

**L00822** 1-Cyclohexyl-3-(2-morpholinoethyl)-carbodiimide methyl-p-toluenesulfonate

**A10973** N,N'-Dicyclohexylcarbodiimide [DCC]

**A19292** N,N'-Diisopropylcarbodiimide

**B25057** 1-(3-Dimethylaminopropyl)-3-ethylcarbodiimide

**A10807** 1-(3-Dimethylaminopropyl)-3-ethylcarbodiimide hydrochloride [EDCI]

**A10962** 1-(3-Dimethylaminopropyl)-3-ethylcarbodiimide methiodide

**L19462** N-Isopropylcarbodiimide on Merrifield resin

### Active ester reagents

**A10802** Acetone oxime

**L00290** Bis(4-nitrophenyl) carbonate

**B24153** 2,5-Diphenyl-4-hydroxy-3-oxo-2,3-dihydrothiophene 1,1-dioxide [HOTDO]

**L13513** 4,6-Diphenylthieno[3,4-d]-1,3-dioxol-2-one 5,5-dioxide [TDO activated carbonate]

**A12892** 1-Hydroxybenzotriazole hydrate [HOBt]

**A13205** endo-N-Hydroxy-5-norbornene-2,3-dicarboximide [HONB]

**A13862** N-Hydroxyphthalimide

**A14522** 2-Hydroxypyridine

**A10312** N-Hydroxysuccinimide

**A15312** 3-Nitrophenol

**A14376** 4-Nitrophenol

**L00359** 4-Nitrophenyl trifluoroacetate

**A15574** Pentafluorophenol [PFP]

**B25671** 2,4,5-Trichlorophenol

**A10788** 2,2,2-Trifluoroethanol

## Alfa Aesar Worldwide Sales Offices

### NORTH AMERICA

Tel: 1-800-343-0660 or  
1-978-521-6300

Fax: 1-800-322-6350  
Email: info@alfa.com

### GERMANY

Tel: 00800 4566 4566 or  
+49-721-84007-280

Fax: 00800 4577 4577 or  
Email: Eurosales@alfa.com

### UNITED KINGDOM

Tel: 0800-801812 or  
+44 (0)1524-850506

Fax: +44 (0)1524-850608  
Email: UKsales@alfa.com

### FRANCE

Tel: 0800 03 51 47 or  
+33 03 88 62 26 90

Fax: +33 03 88 62 26 81  
Email: frventes@alfa.com

### INDIA

Tel: +91 (0)44 3251 9881  
Fax: +91 (0)44 2814 4154

Email: IndiaSales@alfa.com

### CHINA

Tel: (010) 51660581  
Fax: (010) 51987480

Email: sales-china@alfa.com

# Peptide Reagents Con't

## Pyridinium salts

**A12820** 2-Chloro-1-methylpyridinium iodide

**L11088** 2-Fluoro-1-methylpyridinium p-toluenesulfonate

## Phosponium salts

**A16140** Benzotriazol-1-yloxytris-(dimethylamino) phosphonium hexafluorophosphate [BOP]

**L19384** Bromotri(pyrrolidino)phosphonium hexafluorophosphate [PyBroP]

**B25251** Benzotriazol-1-yloxytris(pyrrolidino) phosphonium hexafluorophosphate [PyBOP]

## Uronium salts

**B23597** O-(1H-Benzotriazol-1-yl)-N,N,N',N'-tetramethyluronium hexafluorophosphate [HBTU]

**L18496** O-[(Ethoxycarbonyl)cyanomethylene-amino]-N,N,N',N'-tetramethyluronium tetrafluoroborate [TOTU]

**L13470** O-(1H-Benzotriazol-1-yl)-N,N,N',N'-tetramethyluronium tetrafluoroborate [TBTU]

**L13538** 2-Succinimido-1,1,3,3,-tetramethyluronium tetrafluoroborate [TSTU]

**L19494** O-(1,2-Dihydro-2-oxo-1-pyridyl)-N,N,N',N'-tetramethyluronium tetrafluoroborate [TPTU]

## Miscellaneous coupling reagents

**L14107** Diethyl cyanophosphonate

**A15087** Ethyl diphenylphosphinite

**A12124** Diphenylphosphonic azide

**L14159** 2,4,6-Trichlorobenzoyl chloride

**A11118** 2,2'-Dipyridyl disulfide

## BASES

Most coupling reactions require the presence of a tertiary amine base such as triethylamine. Milder bases, especially 4-methylmorpholine, are less likely to promote racemization and other side reactions. Non-nucleophilic bases such as N-ethyldiisopropylamine (Hünig's Base) are also widely used. The 4-dialkylaminopyridines are frequently used in substoichiometric amounts, in combination with another amine, as hyper-nucleophilic catalysts in mixed anhydride or carbodiimide coupling reactions. Tertiary amine bases are often also required in the introduction of protecting groups.

## Tertiary amine bases

**A11058** 2,4,6-Collidine

**A12158** 4-Methylmorpholine

**L14143** 2,6-Di-tert-butyl-4-methylpyridine

**L03398** 1-Methylpiperidine

**L05265** 1-Diethylamino-2-propanol

**L19372** Morpholine, polymer supported [Methylmorpholine on polystyrene]

**A11081** N-Ethyldiisopropylamine

**A11905** 4-Ethylmorpholine

**A12314** 1,1,3,3-Tetramethylguanidine

**A16294** 1-Ethylpiperidine

**B23797** Tribenzylamine

**A10478** 2,6-Lutidine

**A12646** Triethylamine

## Acylation catalysts

**A13016** 4-Dimethylaminopyridine [DMAP]

**A12575** 1-Methylimidazole

**B23989** 4-(1-Pyrrolidino)pyridine

**A11597** 1,2,4-Triazole

## Alfa Aesar Worldwide Sales Offices

### NORTH AMERICA

Tel: 1-800-343-0660 or  
1-978-521-6300  
Fax: 1-800-322-6350  
Email: info@alfa.com

### GERMANY

Tel: 00800 4566 4566 or  
+49-721-84007-280  
Fax: 00800 4577 4577 or  
Email: Eurosales@alfa.com

### UNITED KINGDOM

Tel: 0800-801812 or  
+44 (0)1524-850506  
Fax: +44 (0)1524-850608  
Email: UKSales@alfa.com

### FRANCE

Tel: 0800 03 51 47 or  
+33 03 88 62 26 90  
Fax: +33 03 88 62 26 81  
Email: frventes@alfa.com

### INDIA

Tel: +91 (0)44 3251 9881  
Fax: +91 (0)44 2814 4154  
Email: IndiaSales@alfa.com

### CHINA

Tel: (010) 51660581  
Fax: (010) 51987480  
Email: sales-china@alfa.com

# Peptide Reagents Con't

## MISCELLANEOUS REAGENTS AND SOLVENTS

The Alfa Aesar product range includes many other items useful in peptide chemistry. This section lists a limited number of examples of materials not covered by earlier sections of this Appendix. The Catalogue also contains a comprehensive listing of natural and unnatural amino acids and their protected derivatives.

### Deprotecting reagents

The Catalog text entries for the protecting reagents listed in this Appendix give some details of methods of cleavage. Some of the most frequently-used reagents are listed below. For further information, see the general references.<sup>1-15</sup>

- |   |   |
|---|---|
| <b>L11577</b> 4-(Aminomethyl)piperidine               | <b>42578</b> Palladium hydroxide on carbon<br>[Pearlman's Catalyst] |
| <b>A15275</b> Boron trifluoride diethyl ether complex | <b>A12442</b> Piperidine  |
| <b>A14740</b> Dichloroacetic acid                     | <b>L13303</b> Tetra-n-butylammonium fluoride trihydrate             |
| <b>A14005</b> Hydrazine monohydrate                   | <b>L06374</b> Trifluoroacetic acid, 99%                             |
| <b>A14475</b> Hydrobromic acid, 45% in acetic acid    | <b>A12198</b> Trifluoroacetic acid, 99%                             |
| <b>L17117</b> Hydrogen fluoride pyridine complex      | <b>A14365</b> Trifluoroacetic acid, biochem. grade, 9.5+%           |
| <b>A10355</b> Morpholine                              | <b>A10173</b> Trifluoromethanesulfonic acid                         |
| <b>A13565</b> Methanesulfonic acid                    | <b>B21789</b> Tris(2-aminoethyl)amine                               |
| <b>A12623</b> Palladium, 5% on carbon                 | <b>L19373</b> Tris(2-aminoethyl)amine, polymer-supported            |
| <b>A12012</b> Palladium, 10% on carbon                |   |

### Cation scavengers

Cleavage of protecting groups under acidic conditions often liberates cations (for example tert-butyl, benzyl, or trityl ions) which can undergo unwanted side reactions with the peptide molecule. The following compounds are among the most useful cation scavengers for these reactive species:

- |                                  |                                  |
|----------------------------------|----------------------------------|
| <b>A12997</b> Anisole            | <b>A14846</b> Thioanisole        |
| <b>L12865</b> Ethanedithiol      | <b>A10320</b> Triethylsilane     |
| <b>L04163</b> Pentamethylbenzene | <b>L09585</b> Triisopropylsilane |

### Solvents

A wide variety of solvents can be used in peptide synthesis. Those listed are among the most useful. Lack of solubility often dictates the use of relatively polar, including dipolar aprotic, solvents. Solvent mixtures can also be used:

- |                                     |   |
|-------------------------------------|---|
| <b>A19862</b> Acetonitrile          | <b>A12747</b> 1,1,1,3,3,3-Hexafluoro-2-propanol |
| <b>L13089</b> Dichloromethane       | <b>A12260</b> 1-Methyl-2-pyrrolidinone          |
| <b>A10924</b> N,N-Dimethylacetamide | <b>L13304</b> Tetrahydrofuran                   |
| <b>A13547</b> N,N-Dimethylformamide | <b>A10788</b> Trifluoroethanol                  |
| <b>A13280</b> Dimethyl sulfoxide    |   |

## Alfa Aesar Worldwide Sales Offices

### NORTH AMERICA

Tel: 1-800-343-0660 or  
1-978-521-6300  
Fax: 1-800-322-6350  
Email: info@alfa.com

### GERMANY

Tel: 00800 4566 4566 or  
+49-721-84007-280  
Fax: 00800 4577 4577 or  
Email: Eurosales@alfa.com

### UNITED KINGDOM

Tel: 0800-801812 or  
+44 (0)1524-850506  
Fax: +44 (0)1524-850608  
Email: UKsales@alfa.com

### FRANCE

Tel: 0800 03 51 47 or  
+33 03 88 62 26 90  
Fax: +33 03 88 62 26 81  
Email: frventes@alfa.com

### INDIA

Tel: +91 (0)44 3251 9881  
Fax: +91 (0)44 2814 4154  
Email: IndiaSales@alfa.com

### CHINA

Tel: (010) 51660581  
Fax: (010) 51987480  
Email: sales-china@alfa.com

# Peptide Reagents Con't



## Miscellaneous reagents and materials

- |   |   |
|---|---|
| <b>L19593</b> N-Acryloylsarcosine methyl ester resin<br>[Sheppard acrylamide resin] | <b>B22183</b> 4-(Hydroxymethyl)benzoic acid       |
| <b>L19600</b> BT-Core resin   | <b>L15721</b> 4-(Hydroxymethyl)phenoxyacetic acid |
| <b>L19602</b> 2-Chlorotrityl alcohol on polystyrene                                 | <b>L17027</b> Merrifield Resin, 1% crosslinked    |
| <b>12457</b> Copper(II) chloride  | <b>A16087</b> Merrifield Resin, 2% crosslinked    |
| <b>A10138</b> 1,4-Dithioerythritol  | <b>L17028</b> Wang Resin, 1% crosslinked          |
| <b>A15797</b> 1,4-Dithio-DL-threitol  | <b>L19369</b> Wang Resin, 2% crosslinked          |

## References

- (1) P. D. Bailey, *An Introduction to Peptide Chemistry*, Wiley, Chichester (1990).
- (2) J. Jones, *The Chemical Synthesis of Peptides*, OUP, Oxford (1991).
- (3) M. Bodansky, *Peptide Chemistry: a Practical Textbook*, Springer-Verlag, Berlin (1993).
- (4) M. Bodansky, *Principles of Peptide Synthesis*, 2nd ed., Springer-Verlag, Berlin (1993).
- (5) M. Bodansky, A. Bodansky, *Practice of Peptide Synthesis*, 2nd ed., Springer-Verlag, Berlin (1994).
- (6) *Peptides: Synthesis, Structures and Applications*, B. Gutte, Ed., Academic Press, San Diego (1995).
- (7) E. Atherton, R. C. Sheppard, *Solid Phase Peptide Synthesis*, IRL Press, Oxford (1989).
- (8) 'Side reactions in peptide synthesis', M. Bodansky, J. Martinez, *Synthesis*, 333 (1981).
- (9) 'Convergent solid-phase peptide synthesis', P. Lloyd-Williams *et al*, *Tetrahedron*, **49**, 11065 (1993).
- (10) 'Synthesis of peptides with mixed anhydrides', N. F. Albertson, *Org. React.*, **12**, 157 (1962).
- (11) 'The azide method in peptide synthesis', Y. S. Klausner, M. Bodansky, *Synthesis*, 549 (1974).
- (12) 'Peptide synthesis via amino acid halides', L. A. Carpino *et al*, *Acc. Chem. Res.*, **29**, 268 (1996).
- (13) T. W. Greene, P. G. M. Wuts, *Protective Groups in Organic Synthesis*, 3rd ed., Wiley, N.Y. (1999).
- (14) P. J. Kocienski, *Protecting Groups*, Thieme, Stuttgart (1994).
- (15) 'Protecting group strategies in organic synthesis', M. Schelhaas, H. Waldmann, *Angew. Chem. Int. Ed.*, **35**, 2057 (1996).
- (16) 'Bibliographical and critical study of the protection of the thiol function in peptide synthesis', F. Cavelier, J. Daunis, R. Jacquier, *Bull. Soc. Chim. Fr.*, 210 (1990).
- (17) J. C. Sheehan, G. P. Hess, *J. Am. Chem. Soc.*, **77**, 1067 (1955); J. C. Sheehan, P. A. Cruickshank, G. L. Boshart, *J. Org. Chem.*, **26**, 2525 (1961).
- (18) T. Curtius, *Ber.*, **35**, 3226 (1902).
- (19) T. Shiori, K. Ninomiya, S. Yamada, *J. Am. Chem. Soc.*, **94**, 6203 (1972).

## Alfa Aesar Worldwide Sales Offices

### NORTH AMERICA

Tel: 1-800-343-0660 or  
1-978-521-6300  
Fax: 1-800-322-6350  
Email: info@alfa.com

### GERMANY

Tel: 00800 4566 4566 or  
+49-721-84007-280  
Fax: 00800 4577 4577 or  
Email: Eurosales@alfa.com

### UNITED KINGDOM

Tel: 0800-801812 or  
+44 (0)1524-850506  
Fax: +44 (0)1524-850608  
Email: UKsales@alfa.com

### FRANCE

Tel: 0800 03 51 47 or  
+33 03 88 62 26 90  
Fax: +33 03 88 62 26 81  
Email: frventes@alfa.com

### INDIA

Tel: +91 (0)44 3251 9881  
Fax: +91 (0)44 2814 4154  
Email: IndiaSales@alfa.com

### CHINA

Tel: (010) 51660581  
Fax: (010) 51987480  
Email: sales-china@alfa.com