


# Methanesulfinic acid sodium salt: the high purity alternative




The methyl sulfone group is present in a number of pharmacologically active molecules. Traditional methods of synthesis often involve nucleophilic attack with the malodorous methylmercaptide, followed by oxidation. Alfa Aesar offers Methanesulfinic acid sodium salt, 97%, which can react with electrophilic carbon atoms to form methyl sulfones in a single step.



#### References:

- 1 D. A. Conlon, *et al*, *Org. Process Res. Dev.*, 2006, **10**, 36.
- 2 W. J. Hennen, C.-H. Wong, *J. Org. Chem.*, 1989, **54**, 4692.
- 3 K. R. Campos, M. Jourmet, S. Lee, E. J. J. Grabowski, R. D. Tillyer, *J. Org. Chem.*, 2004, **70**, 268.
- 4 C. Beaulieu, D. Guay, Z. Wang, D. A. Evans, *Tetrahedron Lett*, 2004, **45**, 3233.



For information on these and over 2000 other new products,  
look in the new "All in One" catalog.

[www.alfa.com](http://www.alfa.com)

#### Alfa Aesar Worldwide Sales Offices

##### NORTH AMERICA

Tel: 1-800-343-0660 or  
1-978-521-6300  
Fax: 1-800-322-6350  
Email: [info@alfa.com](mailto:info@alfa.com)

##### GERMANY

Tel: 00800 4566 4566 or  
+49-721-84007-280  
Fax: 00800 4577 4577 or  
Email: [Eurosales@alfa.com](mailto:Eurosales@alfa.com)

##### UNITED KINGDOM

Tel: 0800-801812 or  
+44 (0)1524-850506  
Fax: +44 (0)1524-850608  
Email: [UKsales@alfa.com](mailto:UKsales@alfa.com)

##### FRANCE

Tel: 0800 03 51 47 or  
+33 03 88 62 26 90  
Fax: +33 03 88 62 26 81  
Email: [frventes@alfa.com](mailto:frventes@alfa.com)

##### INDIA

Tel: +91 (0)44 3251 9881  
Fax: +91 (0)44 2814 4154  
Email: [IndiaSales@alfa.com](mailto:IndiaSales@alfa.com)

##### CHINA

Tel: (010) 51660581  
Fax: (010) 51987480  
Email: [sales-china@alfa.com](mailto:sales-china@alfa.com)